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- 1、Xin mei Duan,Xu Yang,Nianlian Mou,et al.Reactive Oxygen Species-Triggered Carbon Dots-Based Biomimetic Nanotheranostic Agents for Atherosclerosis Management.ADVANCED FUNCTIONAL MATERIALS.July 2024.10.1002/adfm.202405629.(IF 18.5)
 - 2、Mengyue Wang,Xun Zhou,Siyi Zhou,et al.Mechanical force drives the initial mesenchymal-epithelial interaction during skin organoid development.Theranostics.May 2023.10.7150/thno.83217.(IF 12.4)
 - 3、Siyi Zhou,Zeming Li,Xinzhong Li,et al.Crosstalk between endothelial cells and dermal papilla entails hair regeneration and angiogenesis during aging.Journal of Advanced Research.May 2024.10.1016/j.jare.2024.05.006.(IF 11.4)
 - 4、Zhu xian Wang,Yi Hu,Yaqi Xue,et al.4'-OH as the Action Site of Lipids and MRP1 for Enhanced Transdermal Delivery of Flavonoids.ACS Applied Materials & Interfaces.March 2023.10.1021/acsami.2c18086.(IF 9.5)
 - 5、Jialu Luan,Liang Wen,Yehua Bao,et al.Systemic toxicity of biodegradable polyglycolic acid microplastics on the gut-liver-brain axis in zebrafish.SCIENCE OF THE TOTAL ENVIRONMENT.October 2024.10.1016/j.scitotenv.2024.176898.(IF 8.2)
 - 6、Tian Zhu,Chunxia Yao,Shuang Hong,et al.Multi-omics reveal toxicity mechanisms underpinning nanoplastic in redclaw crayfish (Cherax quadricarinatus).SCIENCE OF THE TOTAL ENVIRONMENT.August 2024.10.1016/j.scitotenv.2024.175680.(IF 8.2)
 - 7、Mengying Yang,Jinyu Qiu,Xin Zhao,et al.6-benzylaminopurine exposure induced development toxicity and behaviour alteration in zebrafish (Danio rerio).ENVIRONMENTAL POLLUTION.March 2021.10.1016/j.envpol.2021.116887.(IF 8.071)
 - 8、Qing Tang,Xi Wang,Qichun Zhou,et al.Fuzheng Kang-Ai inhibits NSCLC cell proliferation via regulating hsa_circ_0048091/hsa-miR-378g/ARRDC3 pathway.PHYTOMEDICINE.April 2023.10.1016/j.phymed.2023.154819.(IF 7.9)
 - 9、Shuhui Zhang,Ping Mi,Jialu Luan,et al.Flourene-9-bisphenol acts on the gut-brain axis by regulating oxytocin signaling to disturb social behaviors in zebrafish.ENVIRONMENTAL RESEARCH.May 2024.10.1016/j.envres.2024.119169.(IF 7.7)
 - 10、Xu Zhang,Jian Li,Jin Chen,et al.Enhanced bone regeneration via PHA scaffolds coated with polydopamine-captured BMP2.Journal of Materials Chemistry B.July 2022.10.1039/D2TB01122K.(IF 7.571)
 - 11、Lili Fan,Xueying Zhang,Yihai Huang,et al.Homoplantagin attenuates high glucose-induced vascular endothelial cell apoptosis through promoting autophagy via the AMPK/TFEB pathway.PHYTOTHERAPY RESEARCH.March 2023.10.1002/ptr.7797.(IF 7.2)
 - 12、Shuhui Zhang,Zhenhua Fu,Yixin Xu,et al.The masculinization steroid milieu caused by fluorene-9-bisphenol disrupts sex-typical courtship behavior in female zebrafish (Danio rerio).ECOTOXICOLOGY AND ENVIRONMENTAL SAFETY.October 2022.10.1016/j.ecoenv.2022.114174.(IF 7.129)
 - 13、Yixin Xu,Shuhui Zhang,Shaozhi Zhang,et al.Exposure of zebrafish embryos to sodium propionate disrupts circadian behavior and glucose metabolism-related development. ECOTOXICOLOGY AND ENVIRONMENTAL SAFETY.June 2022.10.1016/j.ecoenv.2022.113791.(IF 7.129)
 - 14、Shuhui Zhang,Xiang Zuo,Jialu Luan,et al.The deleterious effects and potential therapeutic strategy of fluorene-9-bisphenol on circadian activity and liver diseases in zebrafish and mice.Journal of Environmental Sciences.September 2023.10.1016/j.jes.2023.09.014.(IF 6.9)
 - 15、Jialu Luan,Shuhui Zhang,Yixin Xu,et al.Effects of microplastic exposure on the early developmental period and circadian rhythm of zebrafish (Danio rerio): A comparative study of polylactic acid and polyglycolic acid. ECOTOXICOLOGY AND ENVIRONMENTAL SAFETY.May 2023.10.1016/j.ecoenv.2023.114994.(IF 6.8)
 - 16、Shuhui Zhang,Yixin Xu,Shaozhi Zhang,et al.Flourene-9-bisphenol exposure decreases locomotor activity and induces lipid-metabolism disorders by impairing fatty acid oxidation in zebrafish.LIFE SCIENCES.February 2022.10.1016/j.lfs.2022.120379.(IF 6.78)
 - 17、Mengying Yang,Jiaxing Huang,Shuhui Zhang,et al.Melatonin mitigated circadian disruption and cardiovascular toxicity caused by 6-benzylaminopurine exposure in zebrafish. ECOTOXICOLOGY AND ENVIRONMENTAL SAFETY.July 2021.10.1016/j.ecoenv.2021.112555.(IF 6.291)
 - 18、Mengying Yang,Jialu Luan,Yixin Xu,et al.Cardiotoxicity of Zebrafish Induced by 6-Benzylaminopurine Exposure and Its Mechanism.INTERNATIONAL JOURNAL OF MOLECULAR SCIENCES.July 2022.10.3390/ijms23158438.(IF 6.208)
 - 19、Xinping Guo,Ziang Qian,Su Jiang,et al.Assessing the ecotoxicity of florfenicol exposure at environmental levels: A case study of histology, apoptosis and microbiota in hepatopancreas of Eriocheir sinensis. ECOT
-

-
- OXICOLOGY AND ENVIRONMENTAL SAFETY. February 2024. 10.1016/j.ecoenv.2024.116092. (IF 6.2)
- 20、Wan Huan, You Ting, et al. circ_0003204 Regulates Cell Growth, Oxidative Stress, and Inflammation in ox-LDL-Induced Vascular Endothelial Cells via Regulating miR-942-5p/HDAC9 Axis. *Frontiers in Cardiovascular Medicine*. April 2021. 10.3389/fcvm.2021.646832. (IF 6.05)
 - 21、Kai hui Xie, Zunqiang Yan, Qiaoli Yang, et al. Inc001776 Affects CPB2 Toxin-Induced Excessive Injury of Porcine Intestinal Epithelial Cells via Activating JNK/NF- κ B Pathway through ssc-let-7i-5p/IL-6 Axis. *Cells*. March 2023. 10.3390/cells12071036. (IF 6)
 - 22、Liu Mengqian, Wei Xuerong, Zheng Zijun, et al. AMPK activation eliminates senescent cells in diabetic wound by inducing NCOA4 mediated ferritinophagy. *MOLECULAR MEDICINE*. May 2024. 10.1186/s10020-024-00825-8. (IF 6)
 - 23、Pin g Mi, Qian Gao, Ze-Yang Feng, et al. Melatonin attenuates 17 β -trenbolone induced insomnia-like phenotype and movement deficiency in zebrafish. *CHEMOSPHERE*. April 2020. 10.1016/j.chemosphere.2020.126762. (IF 5.778)
 - 24、Gui ying Sun, Hua Ye, Huijuan Liu, et al. ZPR1 is an immunodiagnostic biomarker and promotes tumor progression in esophageal squamous cell carcinoma. *CANCER SCIENCE*. November 2023. 10.1111/cas.16004. (IF 5.7)
 - 25、Wan g Ting, Li Chaopeng, Shi Min, et al. Circular RNA circZNF532 facilitates angiogenesis and inflammation in diabetic retinopathy via regulating miR-1243/CARM1 axis. *Diabetology & Metabolic Syndrome*. -01 2022. 10.1186/s13098-022-00787-z. (IF 5.395)
 - 26、Hua ng Mei, Guo Tian, Meng Yan, et al. Comprehensive analysis of the prognosis and immune effect of the oncogenic protein Four Jointed Box 1. *Frontiers in Oncology*. May 2023. 10.3389/fonc.2023.1170482. (IF 4.7)
 - 27、Lia ng Yi, Lei Ouyang, Shuang Wang, et al. Long noncoding RNA PTPRG-AS1 acts as a microRNA-194-3p sponge to regulate radiosensitivity and metastasis of nasopharyngeal carcinoma cells via PRC1. *JOURNAL OF CELLULAR PHYSIOLOGY*. April 2019. 10.1002/jcp.28547. (IF 4.522)
 - 28、Xin Shen, Hai Shi, Xinli Chen, et al. Esculetin Alleviates Inflammation, Oxidative Stress and Apoptosis in Intestinal Ischemia/Reperfusion Injury via Targeting SIRT3/AMPK/mTOR Signaling and Regulating Autophagy. *Journal of Inflammation Research*. August 2023. 10.2147/JIR.S413941. (IF 4.5)
 - 29、Na Wang, Yunfeng Guo, Liqin Song, et al. Circular RNA intraflagellar transport 80 facilitates endometrial cancer progression through modulating miR-545-3p/FAM98A signaling. *Journal of Gynecologic Oncology*. October 2021. 10.3802/jgo.2022.33.e2. (IF 4.401)
 - 30、Qia n Tian, Huo Bengang, Deng Xiaorong, et al. Decreased TAX1BP1 participates in systemic lupus erythematosus by regulating monocyte/macrophage function. *INTERNATIONAL IMMUNOLOGY*. July 2023. 10.1093/intimm/dxad027. (IF 4.4)
 - 31、Yu Wang, Xiangyun Chen, Qiduan Chen, et al. Expression of human dCTP pyrophosphatase 1 (DCTPP1) and its association with cisplatin resistance characteristics in ovarian cancer. *JOURNAL OF CELLULAR AND MOLECULAR MEDICINE*. April 2024. 10.1111/jcmm.18371. (IF 4.3)
 - 32、Ma Junfu, Meng Qingliang, Zhan Junping, et al. Paeoniflorin Suppresses Rheumatoid Arthritis Development via Modulating the Circ-FAM120A/miR-671-5p/MDM4 Axis. *INFLAMMATION*. August 2021. 10.1007/s10753-021-01504-0. (IF 4.092)
 - 33、Zhe ng Juyan, Zhang Wei, Kang PeiPei, et al. Midazolam Ameliorates Impairment of the Blood-Brain Barrier (BBB) Against LPS. *NEUROTOXICITY RESEARCH*. April 2022. 10.1007/s12640-022-00508-4. (IF 3.978)
 - 34、Wu Yao, Bai Yu, Lu Yan, et al. Transcriptome sequencing and network pharmacology-based approach to reveal the effect and mechanism of Ji Chuan Jian against Parkinson's disease. *BMC Complementary Medicine and Therapies*. June 2023. 10.1186/s12906-023-03999-6. (IF 3.9)
 - 35、Lia ng Wen, Xiaoting Man, Jialu Luan, et al. Early-life exposure to five biodegradable plastics impairs eye development and visually-mediated behavior through disturbing hypothalamus-pituitary-thyroid (HPT) axis in zebrafish larvae. *COMPARATIVE BIOCHEMISTRY AND PHYSIOLOGY C-TOXICOLOGY & PHARMACOLOGY*. July 2024. 10.1016/j.cbpc.2024.109981. (IF 3.9)
 - 36、Pen g Chu, Sheng Luo, Haozhe Wang, et al. Interactive effects of temperature and salinity on the apoptosis, antioxidant enzymes, and MAPK signaling pathway of juvenile pufferfish (*Takifugu fasciatus*). *Aquaculture Reports*. February 2023. 10.1016/j.aqrep.2023.101483. (IF 3.7)
 - 37、Cha ng Liu, Jin-Ou Wang, Wen-Yang Zhou, et al. Long non-coding RNA LINC01207 silencing suppresses AG R2 expression to facilitate autophagy and apoptosis of pancreatic cancer cells by sponging miR-143-5p. *MOLECULAR AND CELLULAR ENDOCRINOLOGY*. April 2019. 10.1016/j.mce.2019.04.004. (IF 3.693)
 - 38、Wan g Shiliang, Li Wenlu, Yang Long, et al. CircPVT1 facilitates the progression of oral squamous cell carcinoma by regulating miR-143-3p/SLC7A11 axis through MAPK signaling pathway. *FUNCTIONAL & INTEGR*
-

-
- ACTIVE GENOMICS.May 2022.10.1007/s10142-022-00865-5.(IF 3.674)
- 39、Lil i Zhu,Dongwu Shi,Jianghong Cao,et al.LncRNA CASC2 Alleviates Sepsis-induced Acute Lung Injury by Regulating the miR-152-3p/PDK4 Axis.IMMUNOLOGICAL INVESTIGATIONS.June 2021.10.1080/08820139.2021.1928693.(IF 3.657)
- 40、Wan g Xiaoyu,Liu Lili,Zhang Lei,et al.Circ_0057583 facilitates brain microvascular endothelial cell injury through modulating miR-204-5p/NR4A1 axis.METABOLIC BRAIN DISEASE.November 2021.10.1007/s11011-021-00866-6.(IF 3.584)
- 41、Yan g Ling,Wang Lijuan,Wang Juan,et al.Long non-coding RNA Gm11974 aggravates oxygen-glucose deprivation-induced injury via miR-122-5p/SEMA3A axis in ischaemic stroke.METABOLIC BRAIN DISEASE.August 2021.10.1007/s11011-021-00792-7.(IF 3.584)
- 42、Kai hui Xie,Qiaoli Yang,Zunqiang Yan,et al.Identification of a Novel lncRNA LNC_001186 and Its Effects on CPB2 Toxin-Induced Apoptosis of IPEC-J2 Cells.Genes.May 2023.10.3390/genes14051047.(IF 3.5)
- 43、Jie Gao,Chengzhi Ding,Junhui Zhou,et al.Propofol suppresses lung cancer tumorigenesis by modulating the circ-ERBB2/miR-7-5p/FOXO1 axis.Thoracic Cancer.-01 2021.10.1111/1759-7714.13856.(IF 3.5)
- 44、Min gming Gu,Kun Liu,Hui Xiong,et al.MiR-130a-3p inhibits endothelial inflammation by regulating the expression of MAPK8 in endothelial cells.Heliyon.-01 2024.10.1016/j.heliyon.2024.e24541.(IF 3.4)
- 45、Chen Feng,Jiang Ruilai,et al.Circular RNA circ_0002984 Facilitates the Proliferation and Migration of Ox-LDL-Induced Vascular Smooth Muscle Cells via the Let-7a-5p/KLF5 Pathway.Cardiovascular Toxicology.August 2024.10.1007/s12012-024-09911-z.(IF 3.4)
- 46、Jiu mei Cai,Xiaoping Wang,Wei Liao,et al.Silencing lncRNA 93358 Inhibits the Apoptosis of Myocardial Cells in Myocardial Infarction Rats by Inducing the Expression of SLC8A1.Biomed Research International.July 2022.10.1155/2022/1138709.(IF 3.246)
- 47、Min jie He,Liang Tu,Ruo Shu,et al.Long Noncoding RNA CASC2 Facilitated Wound Healing through miRNA-155/HIF-1 α in Diabetic Foot Ulcers.Contrast Media & Molecular Imaging.June 2022.10.1155/2022/6291497.(IF 3.009)
- 48、Jia Song,Jia Li,Xiaochen Pei,et al.Identification of cuproptosis-related key genes and pathways in Parkinson's disease via bioinformatics analysis.PLoS One.April 2024.10.1371/journal.pone.0299898.(IF 2.9)
- 49、Jia n-Hua Liu,Chen Li,Liang Cao,et al.Exosomal miR-132-3p from mesenchymal stem cells alleviated LPS-induced acute lung injury by repressing TRAF6.AUTOIMMUNITY.September 2021.10.1080/08916934.2021.1966768.(IF 2.815)
- 50、Lin Zenghai,et al.Circ_0004585 Facilitates Tumorigenesis of Colorectal Cancer Via Modulating the miR-38-3p/ZFX Axis and Activating the MEK/ERK Pathway.Cellular and Molecular Bioengineering.-01 2023.10.1007/s12195-022-00756-6.(IF 2.8)

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