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- 1、 Zhan Xu,Jialiang Fan,Weiguo Tian,et al.Cellulose-Based pH-Responsive Janus Dressing with Unidirectional Moisture Drainage for Exudate Management and Diabetic Wounds Healing.ADVANCED FUNCTIONAL MATERIALS.October 2023.10.1002/adfm.202307449.(IF 19)
- 2、 Weijian Yao,Ying Chen,Zehua Li,et al.Single Cell RNA Sequencing Identifies a Unique Inflammatory Macrophage Subset as a Druggable Target for Alleviating Acute Kidney Injury.Advanced Science.February 2022.10.1002/adv.202103675.(IF 17.521)
- 3、 Yuan Jin,Bowen Liu,Qiuyu Li,et al.PAC1 constrains type II inflammation through promotion of CGRP signaling in ILC2s.JOURNAL OF CLINICAL INVESTIGATION.September 2024.10.1172/JCI180109.(IF 13.3)
- 4、 Zhanming Song,Weijian Yao,Xuekang Wang,et al.The novel potential therapeutic target PSMP/MSMP promotes acute kidney injury via CCR2.MOLECULAR THERAPY.May 2024.10.1016/j.ymthe.2024.05.028.(IF 12.1)
- 5、 Zhan Xu,Yuqian Cui,Weiguo Tian,et al.Surface Wrinkled Microsphere Enhanced Irregular Wound Healing Through Synergistic Hygroscopicity,Reversible Wet-Adhesion and Antibacterial Properties.Small Science.January 2024.10.1002/smsc.202300216.(IF 11.1)
- 6、 Zhiyuan Zou,Le Wang,Zhifei Zhou,et al.Simultaneous incorporation of PTH(1–34) and nano-hydroxyapatite into Chitosan/Alginate Hydrogels for efficient bone regeneration.Bioactive Materials.December 2020.10.1016/j.bioactmat.2020.11.021.(IF 8.724)
- 7、 Runshuai Zhang,Yuanzhen He,Luxia Yao,et al.Metformin chlorination byproducts in drinking water exhibit marked toxicities of a potential health concern.ENVIRONMENT INTERNATIONAL.November 2020.10.1016/j.envint.2020.106244.(IF 7.577)
- 8、 Cheng Zhang,Yuyang Chen,Sili Chen,et al.Occurrence,risk assessment,and in vitro and in vivo toxicity of antibiotics in surface water in China.ECOTOXICOLOGY AND ENVIRONMENTAL SAFETY.March 2023.10.1016/j.ecoenv.2023.114817.(IF 6.8)
- 9、 Meng Zhen,Li Zhe,Zhang Erli,et al.Sevelamer Attenuates Bioprosthetic Heart Valve Calcification.Frontiers in Cardiovascular Medicine.September 2021.10.3389/fcvm.2021.740038.(IF 6.05)
- 10、 Huijuan Cao,Keda Shi,Jing Long,et al.PDGF-BB prevents destructive repair and promotes reparative osteogenesis of steroid-associated osteonecrosis of the femoral head in rabbits.BONE.December 2022.10.1016/j.bone.2022.116645.(IF 4.626)
- 11、 Chun Yang,Jing-Jie Zhang,Xiao-Peng Zhang,et al.Sporamin suppresses growth of xenografted colorectal carcinoma in athymic BALB/c mice by inhibiting liver β -catenin and vascular endothelial growth factor expression.WORLD JOURNAL OF GASTROENTEROLOGY.July 2019.10.3748/wjg.v25.i25.3196.(IF 3.411)
- 12、 Long Aizhuan,Zhao Shuyun,Huang Guanyou,et al.The role of the indoleamine 2,3-dioxygenase gene in preventing ovarian transplant rejection in rats.BIOLOGY OF REPRODUCTION.March 2024.10.1093/biolre/roae036.(IF 3.1)
- 13、 Zhengzheng Xiao,Jing Liang,Qing Deng,et al.Pimozide augments bromocriptine lethality in prolactinoma cells and in a xenograft model via the STAT5/cyclin D1 and STAT5/Bcl-xL signaling pathways.INTERNATIONAL JOURNAL OF MOLECULAR MEDICINE.November 2020.10.3892/ijmm.2020.4784.(IF 3.098)
- 14、 Zhengzheng Xiao,Xiaoli Yang,Zebin Liu,et al.GASC1 promotes glioma progression by enhancing NOTCH1 signaling.Molecular Medicine Reports.February 2021.10.3892/mmr.2021.11949.(IF 2.952)
- 15、 Tang Li,Li Yongfeng,Li Ruiqi,et al.Development and structural characteristics of pseudoosteodentine in the Pacific cutlassfish, Trichiurus lepturus.TISSUE & CELL.May 2022.10.1016/j.tice.2022.101847.(IF 2.586)

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