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- 1、 Hongbin Wei, Zhi Song, Yurong Xie, et al.High temperature inhibits vascular development via the PIF 4-miR166-HB15 module in Arabidopsis.CURRENT BIOLOGY.July 2023.10.1016/j.cub.2023.06.049.(IF 9.2)
 - 2、 Zhang Ke, Zhang Hao, Pan Yanyun, et al.Cell- and noncell-autonomous AUXIN RESPONSE FACTOR3 controls meristem proliferation and phyllotactic patterns.PLANT PHYSIOLOGY.August 2022.10.1093/plphys/kiac370.(IF 8.005)
 - 3、 Bai Feng, Ma Huijin, Cai Yichang, et al.Natural allelic variation in GRAIN SIZE AND WEIGHT 3 of wild rice regulates the grain size and weight.PLANT PHYSIOLOGY.May 2023.10.1093/plphys/kiad320.(IF 7.4)
 - 4、 Sheng Xiaoguang, Yu HuiF ang, Wang Jiansheng, et al.Establishment of a stable, effective and universal genetic transformation technique in the diverse species of Brassica oleracea.Frontiers in Plant Science.October 2022.10.3389/fpls.2022.1021669.(IF 6.627)
 - 5、 Min Zheng, Yu Song, Lingyu Wang, et al.CarH57, a RNA helicase, contributes pepper tolerance to heat stress.PLANT PHYSIOLOGY AND BIOCHEMISTRY.November 2023.10.1016/j.plaphy.2023.108202.(IF 6.5)
 - 6、 Qingzhu Wang, Shikang Lei, Jiawen Yan, et al.UBC6, a ubiquitin-conjugating enzyme, participates in secondary cell wall thickening in the inflorescence stem of Arabidopsis.PLANT PHYSIOLOGY AND BIOCHEMISTRY.October 2023.10.1016/j.plaphy.2023.108152.(IF 6.5)
 - 7、 Fan Yongyi, Chen Hongmei, Wang Belf ang, et al.DWARF AND LESS TILLERS ON CHROMOSOME 3 Promotes Tillering in Rice by Sustaining FLORAL ORGAN NUMBER 1 Expression.PLANT PHYSIOLOGY.July 2024.10.1093/plphys/kiac367.(IF 6.5)
 - 8、 Jiawen Yan, Yu Song, Meng Li, et al.IRR1 contributes to de novo root regeneration from Arabidopsis thaliana leaf explants.PHYSIOLÓGIA PLANTARUM.October 2023.10.1111/ppl.14047.(IF 6.4)
 - 9、 Shuping Zhao, Mengying Tan, Yamei Zhu, et al.Combined analysis of microRNA and mRNA profiles provides insights into the pathogenic resistant mechanisms of the lotus rhizome rot.PHYSIOLÓGIA PLANTARUM.October 2023.10.1111/ppl.14045.(IF 6.4)
 - 10、 Jingjing Mao, Jiaping Yuan, Zhijie Mo, et al.Overexpression of NtCBL5A Leads to Necrotic Lesions by Enhancing Na⁺ Sensitivity of Tobacco Leaves Under Salt Stress.Frontiers in Plant Science.September 2021.10.3389/fpls.2021.740976.(IF 5.754)
 - 11、 Huang Jinghua, Qiao Zhengping, Yu Hang, et al.OsRH52A, a DEAD-Box protein, is required for embryo sac development by regulating functional megaspore specification in rice.JOURNAL OF EXPERIMENTAL BOTANY.April 2024.10.1093/jxb/erae180.(IF 5.6)
 - 12、 Min Zheng, Qingzhu Wang, Shikang Lei, et al.PtoMPO1, a negative mediator, functions in poplar drought tolerance.PLANT PHYSIOLOGY AND BIOCHEMISTRY.September 2022.10.1016/j.plaphy.2022.09.004.(IF 5.437)
 - 13、 Danfei Lou, Xinyue Xing, Yunyu Liang, et al.Dendrobine modulates autophagy to alleviate ox-LDL-induced oxidative stress and senescence in HUVECs.DRUG DEVELOPMENT RESEARCH.April 2022.10.1002/ddr.21937.(IF 5.004)
 - 14、 Qian Jie, Li Meng, Zheng Min, et al.Arabidopsis SSB1, a Mitochondrial Single-Stranded DNA-Binding Protein, is Involved in ABA Response and Mitochondrial RNA Splicing.PLANT AND CELL PHYSIOLOGY.August 2021.10.1093/pcp/pcab097.(IF 4.927)
 - 15、 Zhao Chongchong, Li Qihang, Ge Qi, et al.Lectin Receptor-Like Protein Kinase OsNRFG6 is Required for Embryo Sac Development and Fertilization in Neo-Tetraploid Rice.Rice.June 2024.10.1186/s12284-024-00720-0.(IF 4.8)
 - 16、 Xinghai Yang, Baoxuan Nong, Can Chen, et al.OsNPF3.1, a member of the NRT1/PTR family, increases nitrogen use efficiency and biomass production in rice.Crop Journal.July 2022.10.1016/j.cj.2022.07.001.(IF 4.647)
 - 17、 Yuan Guang, Nong Tongjia, Hunpatin Oluwaseyi Setonji, et al.Genome-wide identification of Shaker K⁺ channel family in Nicotiana tabacum and functional analysis of NtSKOR1B in response to salt stress.Frontiers in Plant Science.April 2024.10.3389/fpls.2024.1378738.(IF 4.1)
 - 18、 Sujuan Shi, Lulu An, Jingjing Mao, et al.The CBL-Interacting Protein Kinase NtCIPK23 Positively Regulates Seed Germination and Early Seedling Development in Tobacco (Nicotiana tabacum L.).Plants-Basel.February 2021.10.3390/plants10020323.(IF 3.935)
 - 19、 Min Zheng, Tao Peng, Tingting Yang, et al.Arabidopsis MHP1, a homologue of yeast Mpo1, is involved in ABA signaling.PLANT SCIENCE.October 2020.10.1016/j.plantsci.2020.110732.(IF 3.591)
 - 20、 Zheng Min, Zhu Chunyan, Yang Tingting, et al.GSM2, a transaldolase, contributes to reactive oxygen
-

-
- species homeostasis in Arabidopsis. PLANT MOLECULAR BIOLOGY. June 2020. 10.1007/s11103-020-01022-x. (IF 3.302)
- 21、Shaoying Wen, Jiayu Li, Ziyuan Hao, et al. Overexpression of the LcCUC2-like gene in Arabidopsis thaliana alters the cotyledon morphology and increases rosette leaf number. PeerJ. February 2022. 10.7717/peerj.12615. (IF 3.061)
- 22、Sujuan Shi, Fangzheng Xu, Yuqian Ge, et al. NH_4^+ Toxicity, Which Is Mainly Determined by the High NH_4^+/K^+ Ratio, Is Alleviated by CIPK23 in Arabidopsis. Plants-Basel. April 2020. 10.3390/plants9040501. (IF 2.762)
- 23、Yongzhu Wang, Xiefeng Ye, Kang Yang, et al. Characterization, expression, and functional analysis of polyamine oxidases and their role in selenium-induced hydrogen peroxide production in Brassica rapa. JOURNAL OF THE SCIENCE OF FOOD AND AGRICULTURE. February 2019. 10.1002/jsfa.9638. (IF 2.422)
- 24、Manqing Peng, Yiyun Qin, Yuan He, et al. Cloning and Functional Analysis of PmMYB45, a Transcription Factor in Pinus massoniana. Forests. May 2024. 10.3390/f15050814. (IF 2.4)

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