

苯甲基磺酰氟溶液(PMSF,100mmol/L)

产品简介:

蛋白酶抑制剂由很多种, 包括 EDTA、E-64、 NaVO_3 、亮抑酶肽、抑肽酶、胰蛋白酶抑制剂等均可有效抑制蛋白的降解, 并维持原有的蛋白间相互作用, 苯甲基磺酰氟即 Phenylmethanesulfonyl fluoride, 简称为 PMSF, 分子式为 $\text{C}_7\text{H}_7\text{FO}_2\text{S}$, 分子量为 174.19, 纯度 > 99%。

PMSF 为常用生化试剂, 采用进口试剂配制, 可抑制丝氨酸蛋白酶如胰蛋白酶(trypsin)和胰凝乳蛋白酶(chymotrypsin), 也可抑制半胱氨酸蛋白酶和乙酰胆碱酯酶, 其有效浓度为 0.1 ~ 1mM。该试剂仅用于科研领域, 不适用于临床诊断或其他用途。

产品组成:

名称 \ 编号	PI0011	PI0011	PI0011	Storage
PMSF Solution(100mmol/L)	1.5ml	5ml	10ml	-20°C
使用说明书	1 份			

操作步骤(仅供参考):

- 1、使用前 PMSF 充分溶解混匀, 根据溶液的量加入 PMSF Solution(100mmol/L), 使其终浓度达到 1mM。

注意事项:

- 1、PMSF 有轻微毒性, 小心操作。
- 2、为了您的安全和健康, 请穿实验服并戴一次性手套操作。
- 3、试剂开封后请尽快使用, 以防影响后续实验效果。

有效期: 12 个月有效; 低温运输, -20°C 保存。

相关产品:

产品编号	产品名称
PI0015	蛋白酶抑制剂混合液(100×PIC)
PS0013	RIPA 裂解液(强)
PT0001	BCA 法蛋白定量试剂盒
DH0006	苏木素伊红(HE)染色液(醇溶)
TO1013	丙二醛(MDA)检测试剂盒(TBA 比色法)

文献引用:

- 1、Chen Jianxiong,Duan Shiyu,Wang Yulu,et al.MYG1 drives glycolysis and colorectal cancer development through nuclear-mitochondrial collaboration.Nature Communications.10.1038/s41467-024-49221-0.June 2024.(IF 14.7)
 - 2、Zhiqiang Xie,Liqing Lu,Zhenghua Wang,et al.S-acylthioalkyl ester (SATE)-based prodrugs of deoxyribose cyclic dinucleotides (dCDNs) as the STING agonist for antitumor immunotherapy.EUROPEAN JOURNAL OF MEDICINAL CHEMISTRY.10.1016/j.ejmech.2022.114796.September 2022.(IF 7.088)
 - 3、Zhenxiu Liu,Lin Chen,Mingyun Chen,et al.Sarmentol H derived from Sedum sarmentosum Bunge directly targets FXR to mitigate cholestasis by recruiting SRC-1.PHYTOMEDICINE.10.1016/j.phymed.2024.155759.May 2024.(IF 6.7)
 - 4、Zhou Rui,Lin Chuman,Cheng Yanzhen,et al.Liraglutide Alleviates Hepatic Steatosis and Liver Injury in T2MD Rats via a GLP-1R Dependent AMPK Pathway.Frontiers in Pharmacology.10.3389/fphar.2020.600175.March 2021.(IF 5.811)
 - 5、Xianning Zhang,Lulu Liu,Xin Liu,et al.Chidamide suppresses adipogenic differentiation of bone marrow derived mesenchymal stem cells via increasing REEP2 expression.iScience.10.1016/j.isci.2023.106221.February 2023.(IF 5.8)
 - 6、Zhiqiang Xie,Yuchen Yang,Zhenghua Wang,et al.Dithioethanol (DTE)-Conjugated Deoxyribose Cyclic Dinucleotide Prodrugs (DTE-dCDNs) as STING Agonist.INTERNATIONAL JOURNAL OF MOLECULAR SCIENCES.10.3390/ijms25010086.December 2023.(IF 5.6)
 - 7、Jie Zhou,Ruhe Yang,Yiwen Sun,et al.HCIO-triggered interventional probe enabled early detection and intervention of atherosclerosis.ANALYST.10.1039/D2AN01374F.November 2022.(IF 5.227)
 - 8、Qiaolei Wang,Rui Lou,Qianwen Yin,et al.A nano-detection system based on a chemical probe for early diagnosis of atherosclerosis in situ.ANALYST.10.1039/D1AN00484K.June 2021.(IF 4.616)
 - 9、Yuqi Sang,Chenxu Niu,Jiaxi Xu,et al.PI4KIII β -mediated Phosphoinositides Metabolism Regulates Function of the VTA Dopaminergic Neurons and Depression-like Behavior.JOURNAL OF NEUROSCIENCE.10.1523/JNEUROSCI.0555-23.2024.January 2024.(IF 4.4)
 - 10、Jiahao Yan,Yuhe Ding,Zhao Peng,et al.Systematic Proteomics Study on the Embryonic Development of Danio rerio.JOURNAL OF PROTEOME RESEARCH.10.1021/acs.jproteome.3c00056.July 2023.(IF 4.4)
 - 11、Hang Yu,Rui Feng,Fenglin Chen,et al.Rapid FRET Assay for the Early Detection of Alpha-Synuclein Aggregation in Parkinson's Disease.ACS Chemical Neuroscience.10.1021/acscchemneuro.3c00617.March 2024.(IF 4.1)
 - 12、Shuangli Yang,Chuman Lin,Xiaoyun Zhuo,et al.Glucagon-like peptide-1 alleviates diabetic kidney disease through activation of autophagy by regulating AMP-activated protein kinase-mammalian target of rapamycin pathway.AMERICAN JOURNAL OF PHYSIOLOGY-ENDOCRINOLOGY AND METABOLISM.10.1152/ajpendo.00195.2019.November 2020.(IF 3.469)
 - 13、Minglan Shi,Meinian Xu,Xiaowen Huang,et al.The effect of autophagy on hemoporphin-mediated photodynamic therapy in human umbilical vein endothelial cells.Photodiagnosis and Photodynamic Therapy.10.1016/j.pdpdt.2024.104196.May 2024.(IF 3.1)
 - 14、Dao Ge Chen,Cheng Wu Zhao,Yan Chun Gong,et al.Study on the Influencing Factors of Hypoglycemic Effect of Folate Targeted Polymersomes Encapsulating Insulin.JOURNAL OF PHARMACEUTICAL SCIENCES.10.1016/j.xphs.2020.11.042.December 2020.(IF 2.997)
 - 15、Fangqin Yu,Lele Zhang,Runsheng Ma,et al.The Antitumour Effect of Prunella vulgaris Extract on Thyroid Cancer Cells In Vitro and In Vivo.Evidence-based Complementary and Alternative Medicine.10.1155/2021/8869323.January 2021.(IF 2.63)
 - 16、Jie Yang,Shui-yi Shao,Wan-qing Chen,et al.Cigarette Smoke Extract Exposure: Effects on the Interactions
-

between Titanium Surface and Osteoblasts.Biomed Research International.10.1155/2019/8759568.April 2019.(IF 2.197)

17、Chen Xueqin,Tan Jingyi,Yang Mengqi,et al.Genistein has the function of alleviating and treating disseminated intravascular coagulation caused by lipopolysaccharide.Journal of Natural Medicines.10.1007/s11418-018-1215-9.May 2018.(IF 1.92)

18、Yushan Fang,Liang-Cai Wu,Kanglong Ma,et al.Paeoniflorin alleviates lipopolysaccharide-induced disseminated intravascular coagulation by inhibiting inflammation and coagulation activation.DRUG DEVELOPMENT RESEARCH.10.1002/ddr.21647.February 2020.(IF 1.902)

注：更多使用本产品的文献请参考产品网页